



# Handpiece Repair Prices September 1, 2016

## HIGH SPEED BEARING OVERHAULS

M.W Quiet Air Standard	79.00
M.W. 8000 (New Turbine with 60-Day Warranty)	99.50
M.W. Tradition Standard	79.00
M.W. Tradition-L & Push Button (XGT)	98.00
M.W. Quiet Air Power Lever & Push Button	98.00
Star- 430 and Vista Series Standard	79.00
Star- 430 and Vista Series Push Button	95.00
Star- (LF)	Will Call with Estimate
Kavo Super Torque (Regular & Push Button) 647B/649/635B-95	95.00
Kavo Contact Air (Regular & Push Button)	95.00
Kavo - 635B / 647B / 649B	95.00
MK-Dent / Viper	95.00
NSK NL 9000S	95.00
NSK Alpha Light	85.00
Japan Standard (New Turbine with 60-Day Warranty)	59.00
Lares	95.00
Topair	95.00
Schein Master	79.00
Adec and W & H	95.00
Viper 360	95.00
Vector & MK Dent	95.00

## SLOW SPEED COMPLETE REBUILDS

M.W. Tru-Torc Single Speed	145.00
M.W. Tru-Torc Dual Speed	149.00
M.W. Shorty Single Speed	145.00
M.W. Shorty Dual Speed	149.00
M.W. RDH	Will Call with Estimate
M.W. Rhino	165.00
M.W. Shorty Nosecone	90.00
Star Titan Motor	145.00
Coastal Little Guy and Champion	145.00
Japan Slow Speeds	139.00
Star Titan Nosecone	95.00
Air Scalers and New Shafts	Will Call with Estimate
Kavo Slow Speed &/or Attachments	Will Call with Estimate
Lynx Motor	139.00
Dental Aire	139.00
Medidenta	139.00
Lares	145.00



# Handpiece Repair Prices September 1, 2016

## INSTALLATION OF NEW TURBINES

Champion	132.00
Star- 430 SWL (Auto Chuck)	159.00
Star- 430 SWL (Standard)	99.00
Kavo- 635-B ~ 647-B ~ 649-B~ 659-B~ 6000-B~ 6500-B	195.00
Kavo- 632-B ~ 642-B ~ 645-B ~ 630-B ~ 640-B~ 646-B	175.00
NSK- Mach Light Phatelus (Single Spray Turbines)	197.00
NSK- NL - 85 (Triple Spray Turbines)	218.00
NSK- Alpha Light (Single Spray Turbines)	170.00
NSK - NL 95	237.00
NSK - NL9000S	195.00
Viper	169.00
Midwest Quiet Air & Tradition Standard Chucking	99.00
Midwest Tradition - XGT Push Button	159.00
Midwest Tradition & Quiet Air Power Lever (We can convert to Push Button-below)	165.00
Midwest Push Button Conversion Kit (Includes Turbine & Cap)	185.00
Midwest Stylist	179.00
Midwest 8000 (90 day warranty)	89.00
Japan Canister Push Button with Cap	119.00
Japan Canister Push Button	95.00
Japan Standard Canister	59.00
Lares Replacement Turbines	Varies - Will Call
Kavo Replacement Caps	60.00
Sable Push Button Replacement Caps	40.00
Star Replacement Caps	25.00
Japan Replacement Caps	7.00

## ADDITIONAL REPAIRS

Rethread Connectors	90.00
Kavo Replace Fiber Optic Rod (Including Labor)	Call for Estimate
Kavo Fiber Optic Clip	50.00
Kavo Resolder Head/Ream Head	45.00
Kavo Super Torque 640-B Replace Spray Ring, Retainer Nut & Outer O-Ring	139.00
Fiber Optic Replacement of M.W. Tradition, M.W. Quiet Air, Kavo's 630, 632, 633	145.00
Fiber Optic Replacement of Star 430, NSK Alpha Light, Healthco	145.00
Fiber Optic Replacement of Kavo 640	Call for Estimate
Fiber Optic Replacement of Kavo 642, 643	Call for Estimate
Clear H2O Lines	25.00
Midwest Quite Air - Ream Head	25.00



# Handpiece Repair Prices September 1, 2016

## WARRANTY PERIOD

TURBINE REPLACEMENTS	
All Turbine Replacement <i>(except for canister turbine which are 2 months)</i>	6 Month Warranty
BEARING OVERHAULS	
M.W. 8000 and M.W. 8000-I (Fiber-Optic)	60 Day Warranty
Star Push Button	90 Day Warranty
Kavo - All Models	90 Day Warranty
Midwest Quiet Air and Tradition	90 Day Warranty
Lares	90 Day Warranty
W & H and Adec	90 Day Warranty
MK Dent	90 Day Warranty
Vector	90 Day Warranty

## COMMENTS

Size of Burrs: Recommended for ALL Push Button, Power Lever and Chucks (.0628 - .0630)

Excessive PSI - Check Mfg. specifications for your Handpiece. This varies with different Handpieces. When exceeding recommended pressure, Bearings will show scorching which voids Warranty.

We are unable to warranty water lines due to the fact that debris from Bio Film in system can cause lines to clog at any time.

Call for pricing of Loose Turbines. We are unable to warranty Loose Turbine; Failure is most often caused by improper cleaning and or installation.